


Instruction Sheet

Sealing and Removing Locks and Links



Aluminum Padlock Seal with Reusable Body.

(No Key is Necessary!)

To seal (left), put link through emergency kit eyelet; push padlock up until it locks.

To remove (center), break link at joint.

To replace (right), pull ends of link out of padlock; insert new link.


- 484101 Safety Set (300 Links 100 Bases)
- 484102 Safety Locks (100 Bases)
- 484103 Links (300 each)

- 484111 Mini Pak Set (100 Links 50 Bases)
- 484112 Mini Pak Locks (50 Bases)
- 484113 Mini Pak Links (100 each)

Safety Control Seals

Once locked, it must be broken to remove. No scissors or tools needed to break seal.

Pull to break!



- 484107 Safety Control Seals (Numbered 100/Pkg)
- 484108 Safety Control Seals (Without Numbers 100/Pkg)



- ✉ 800 Glen Avenue, Moorestown, New Jersey 08057
- ☎ 1.800.257.BEAM (2326) • Local: 1.856.359.2231
- ✉ Fax: 1.856.359.2249 • E-mail: info@omnimedbeam.com

Visit our web page on the internet: www.omnimedbeam.com

Products to Help People...Help People